Use the following instructions carefully. 

G-aenial® Universal Flo

UNIVERSAL LIGHT-CURED RADIOPAQUE FLOWABLE COMPOSITE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Direct restorative for Class I, IV, III, II and V cavities.
- Fissure sealant.
- Selecting hypersensitive areas.
- Repair of (or) direct aesthetic restorations, temporary crown & bridge, defect margins when margins are in enamel
- Blocking of microleaks
- Linear or base
- Rehabilitation of crowns & bridges, inlays and veneers using the indirect technique in combination with GC Gradia components (please refer to its technical IFU)
- Splitting of teeth in combination with fibres such as GC EverStick fibres.
- Adhesive luting of ceramic and composite inlays, onlays and veneers with a thickness (< 2.0 mm) that enables the use of a purely light-curing technique.

CONTRAINDICATIONS

Avoid use of this product in patients with known methacrylate monomer and methacrylate polymer allergies

Operation

- Preparations
 - Hold the syringe upright and remove the cap by turning counterclockwise.
 - Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
 - Place the light protective cover until use.
- Note:
 - Fake care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remnants of material.
- Shade Selection
 - Clean the tooth with pumice and water. Shade selection should be made prior to isolation.
 - Select the appropriate shade from 15 shades of A1, A2, A3, A4, B1, B2, B3, C3, A02, A03, BW, CV, JE, AE, A, B, C, shades as based on the Vita® classical shade guide.

- Cavity Preparation
 - Prepare cavity using standard techniques. Dry by gently blowing with oil air. Rub with cotton swab until completely dry. Prepare tooth with contamination with saliva, blood or sultus fluid.
 - Note:
 - For pulp capping, use calcium hydroxide.

- Bonding treatment
 - For bonding G-aenial Universal Flo to enamel and/or dentine, use a light-cured bonding system such as G-aenial® Bond, G-Premio BOND, G-BOND®. Follow manufacturer's instructions.

- Placement of G-aenial Universal Flo
 - Remove cover from the dispensing tip.
 - Note:
 - Prior to syringe material into the cavity, make sure to check the following:
 - Dispensing tip is securely attached to the syringe.
 - Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the tip pointing upwards until material reaches the opening of the tip (Fig. 2).

- Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to syringe material. Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.
- Note:
 - a. If plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe.
 - b. Minimize exposure to ambient light. Ambient light can shorten the manipulation time.
 - c. After use, immediately remove and dispose the dispensing tip and lightly close the syringe with the cap.

- Clinical Hint
 - In direct contact effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This allows the material to spread from the dispensing tip and provides a smooth surface over the restoration.

- Light Curing
 - Light cure the G-aenial Universal Flo using a light curing unit. Keep light curing as close as possible to the surface. Refer to the following chart for irradiation Time and Effective Depth of Cure.

Irradiation Time and Effective Depth of Cure (ISO 4049)	
Shade	Irradiation time 9 sec. (Plasma arc) (2000mW/cm) 10 sec. (GC G-Light) (1200mW/cm) 20 sec. (HalogenLED) (700 mW/cm)
A1, A2, A3, B1, B2, JE, AE	2.0mm
A3.5, A4, B3, C3, BW, CV, A02, A03	1.5mm

- Note:
 - Material should be placed and light cured in layers. For maximum layer thickness, please consult the following table.
 - Lower light intensity may cause insufficient curing or discoloration of the material.
 - For the indirect technique light cure each layer for 30 seconds with GC LABO LIGHT LV-III II or for 10 seconds with GC STEPLIGHT SL-I. Final cure for 3 minutes using GC LABO LIGHT LV-III

- Finishing and Polishing
 - Finish and polish using standard techniques.

CEMENTATION OF INLAYS, ONLAYS AND VENEERS

- Removal of the temporary restoration
 - After removing the provisional restoration, clean the prepared tooth with brush and slurry of pumice to remove any remaining temporary cement. Rinse the prepared tooth with water spray. Titrat the restoration before the cementation and clean thoroughly.
- Treatment of the restoration
 - Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.
- Apply GC Ceramic Primer I or a similar silane coupling agent to the surface of the restoration according to the manufacturer's instructions. Do not rinse it.
- Treatment of tooth preparation
 - Soak the preparation with 35-37% phosphoric acid, rinse and gently dry according to the manufacturer's instructions.
 - Apply GC G-Premio BOND, GC G-aenial BOND or GC G-BOND on prepared tooth surface, leave undisturbed for 10 seconds and dry thoroughly for 5 seconds using maximum air pressure. Refer to the respective instructions for use.
 - Note:
 - When working with rubber dam is recommended.
 - Make sure to blow dry with a maximum air pressure in order to obtain thin film.
 - In case of using other bonding agents, use only those which have thin film thickness to ensure the good fit of the veneer.
- Light cure the bonding agent for the designated time in the instructions for use.
- Luting procedure
 - Select the suitable shade. Apply G-aenial Universal Flo directly to the bonding surface(s) of the restoration and/or tooth surface. Exposure to ambient light should be avoided during application. Seal the restoration exerting moderate pressure and remove excesses using the following techniques:
 - Keep moderate pressure and remove excesses using a brush. Make sure that excess are thoroughly removed from areas of difficult access (e.g. interproximal areas). Tack-cure the veneer for a few seconds and proceed with the removal of remaining excesses using a probe or another appropriate instrument.
 - Tack cure the excesses by waving the light guide of a curing light over the excesses for a few seconds until they reach a solid rubber consistency and remove excesses using a probe or another appropriate instrument.
 - Note: Reaching a good consistency for excess removal depends on the light-curing unit used. It is important to find the best combination of time, intensity, and distance of the curing light in use.

- Through light-cure each surface of the restoration for 30 seconds. Note: To avoid oxygen inhibition layer at the margins of the restoration, cure the margins with a glycerin gel before final light curing procedure.

- Polish the restoration
 - Polish the restoration margins using appropriate instruments for the polishing of composite resins.

SHADES

A1, A2, A3, A3.5, A4, B1, B2, B3, C3, A02, A03, BW, CV, JE, AE

STORAGE

Store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F) away from high temperatures or direct sunlight. (Shelf life : 3 years from date of manufacture)

PACKAGES

- Syringe 3.4g (20mL), 20 dispensing tips (plastic type), 1 light protective cover
- Dispensing tip package:
 - 30 dispensing tips (plastic type), 2 light protective covers
- 30 dispensing tips (needle type), 2 light protective covers

CAUTION

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take care to avoid ingestion of the material.
- Do not mix with other similar products.
- Dispensing tips are for single use only. To prevent cross contamination between patients, do not reuse the dispensing tip. The dispensing tip cannot be sterilized in an autoclave or chemically.
- Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity.
- Wear protective eye glasses during light cure.
- When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of curing dust.
- Avoid getting material on clothing.
- In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
- Do not use G-aenial Universal Flo in combination with eugenol containing materials as eugenol may hinder G-aenial Universal Flo from setting.
- Do not place the dispensing tip with excessive alcohol as this may cause contamination of the paste.
- Dispose of all wastes according to local regulations.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.georgroup.com/download/> They can also be obtained from your supplier.

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