BRILLIANT FLOW UNIVERSAL COMPOSITE

DE	Gebrauchsinformation	3
EN	Instructions for use	6
ES	Instrucciones de uso	9
PT	Instruções de utilização	11
ZH-CN	使用说明	14
ZH-TW	使用說明	17

Definition

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EN

Definition

Brilliant Flow is a nanofilled, radiopaque flowable composite for special indications. Brilliant Flow is cured in a tooth cavity or on a model with all current light-curing lamps (e.g. Coltolux® LED).

Composition

Methacrylates

Barium glass, silanized

Amorphous silica, hydrophobic

Technical data

Complies with ISO 4049:2000

Average filler particle size: 0.6 µm

Range of particle size: 0.04–2.5 µm

Brilliant Flow:

Percentage by volume of total inorganic filler: 42 %

Percentage by weight of total inorganic filler: 63 %

Indications

- Direct fillings of class V (cervical caries, root erosion, wedge-shaped defects)
- Repairs to anterior teeth (class III, IV)
- Small fillings in the posterior region
- Fillings of minimal cavities
- Repairs of composite and ceramic veneers
- Block out of undercuts
- Adhesive cementing of composite and ceramic restorations, if light penetration is possible.
- Preventive resin restorations

Contra-indications

If any allergies exist to any components within Brilliant Flow. If the

application field cannot be isolated during application and curing of Brilliant Flow. For patients with poor oral hygiene.

Side effects

No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.

Interactions with other agents

Agents containing eugenol and/or oil of cloves may affect the polymerization of Brilliant Flow. The use of zinc-oxide-eugenol cements in combination with Brilliant Flow should therefore be avoided. Discolorations may occur when using cationic mouth rinses as well as plaque indicators or chlorhexidin.

Note:

- Light can lead to premature curing. Recap syringe after each use.
- To prevent run off, pull back the syringe plunger.
- For repairs, the treated surface must be roughened first.

General information

Preliminary preparation

Cleansing the tooth

Cleanse tooth to be treated and the adjacent teeth using a brush and prophylaxis paste which does not contain fluoride.

Selection of shade

Selection of shade should be made before isolation. The shade is ideally determined in broad daylight using the Vita™ Shade Guide. Shades reach their final tone 24 hours after curing.

Isolation of the cavity (rubber dam)

Adequate isolation is absolutely necessary for achieving optimum results. The use of HYGENIC® /Roeko Dental Dams are recommended.

Matrix and inter-dental wedges

When working approximately, apply a thin matrix band. Fix the matrix proximally using interdental wedges.

Adhesive system

Apply adhesive system (e.g. Coltène® One Coat Bond SL or A.R.T. BOND) according to the corresponding instructions for use.

Working time

Brilliant Flow is light sensitive and should not be left too long under intensive light before curing, especially operatory light or sun light.

Polymerization

Brilliant Flow can be cured using all current curing lights (e.g. Coltolux® LED).

Recommended exposure times:

Layer Thickness: 2 mm

> 500 mW/cm²

Brilliant Flow A1/B1	20 s
Brilliant Flow A2/B2	20 s
Brilliant Flow A3/D3	20 s
Brilliant Flow A3.5/B3	20 s
Brilliant Flow A4/M5	20 s
Brilliant Flow Super White	20 s
Brilliant Flow Transparent	20 s

Finishing

For the tooth preparation, appropriate rotary instruments and polishers from the Diatech / Alpen range can be used. Fluoridate all surfaces after completion.

Application of Brilliant Flow for direct restorations

Preparation of the cavity

When preparing the cavity, every effort possible should be made to use a technique that preserves tooth structure (principle of the Adhesive Restoration Technique). Prepare enamel and dentin using 80 µm preparation diamonds and finish using 25 µm finishing diamonds. Bevelled enamel margins are recommended to increase the area of adhesion between tooth and filling material and thereby optimize the marginal seal.

Pulp Protection

Covering the areas close to the pulp with hard-setting calcium hydroxide cements using the spot technique protects the pulp against bacterial infiltration.

Application of Brilliant Flow

Apply the material using the mounted applicator tip directly into the cavity.

Application of Brilliant Flow for indirect restorations

Pretreat contact surfaces of cavity according to the corresponding manufacturer's instructions for use.

Application of Brilliant Flow

Dispense the material using the mounted applicator tip directly into the restoration or cavity preparation. Position restoration in place using gentle pressure.

Removal of excess material

Rough excess material (e.g. brush or spatula). Afterwards, hold restoration with pressure in position and remove additional excess material with a spatula or a suitable instrument.

Polymerization

Recommended curing time: cure each side for 40 s.

Emergency measures

In the case of direct contact with the oral mucosa, rinsing with tap water is sufficient. In case of contact with the eyes, rinse thoroughly with water (10 min), and consult an ophthalmologist showing these instructions for use.

Note

Only supplied to dentists and dental laboratories or upon their instructions. Keep out of the reach of children! Should not be used after expiration date.

Shelf life and marking

The expiration date and the **LOT** number are marked on the packaging containers.

Storage

4–23 °C / 39–73 °F. Avoid exposure to direct sunlight or other heat sources.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

Infection control

For hygienic reasons, the application tips from Brilliant Flow are only designed for single use.

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Coltène/Whaledent AG