**GRADIA® DIRECT LoFlo**

LIGHT-CURED FLOWABLE COMPOSITE

GRADIA® DIRECT LoFlo is a multi-purpose light-cured, radiopaque, fluoride releasing, micro-filled hybrid type composite resin. Thanks to GC's new HDR (High Density Radiopaque Pre-Polymerized Filler) technology, it contributes to achieve excellent abrasion resistance. Also the thicker consistency (lower flowability) is excellent for cavities of class III, V and root surface caries.

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

- Restoration of Class I, II, III, IV, V cavities (particularly for small Class I cavities / shallow Class V cavities / other small cavities).
- Restoration of root surface caries.
- Restorations in deciduous teeth.
- Filling tunnel shaped cavities.
- Sealing hypersensitive areas.
- Liner / base / filling in cavity undercuts.
- Sealant.
- Fixing mobile teeth.
- Additions to composite restorations. (Fig. 1)

**CONTRAINDICATIONS**

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**DIRECTIONS FOR USE**

## 1. Preparations

- Hold the syringe upright and remove the wing cap by turning counterclockwise. Take care not to expose material to direct light from the dental lamp or natural light (Fig. 2).
- Promptly and securely attach the dispensing tip (plastic or needle type) to the syringe by turning clockwise (Fig. 2)(Fig. 3).

Note : Take care not to attach the dispensing tip too tightly. This may damage its screw.

- After attaching the dispensing tip, protect it with the cover until ready for use in order to avoid exposure to light (Fig. 2).

## 2. Cavity Preparation

Prepare cavity using standard techniques.

Note : For pulp capping, use calcium hydroxide.

## 3. Shade Selection

Select shade from 7 shades of A1, A2, A3, A3.5, AO3, CV and BW.

## 4. Use of Light-cured Bonding System

For bonding GRADIA® DIRECT LoFlo to tooth structure, use a light-cured bonding system such as GC UniFil® Bond or G-BOND.

Note : When using light-cured bonding system, follow manufacturer's instructions for use.

## 5. Placement of GRADIA® DIRECT LoFlo

- Remove cover from the dispensing tip on the syringe (Fig. 5).

2) Prior to extrusion of material, gently trial push the syringe plunger outside the mouth to make sure that the dispensing tip is securely attached to the syringe. To remove any air from the dispensing tip, with the tip pointing upwards gently push forward the syringe plunger until material reaches the mouth of the tip (Fig. 6).

Note : If there is air inside the dispensing tip, air bubbles may be formed at the time of injection.

- Place the dispensing tip as close as possible to the cavity, and slowly push the plunger to inject material into it (Fig. 7). Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

Note :

- When attaching the dispensing tip, make sure that no material is sticking to the joint between the tip and the syringe in order to ensure a tight connection.

b. If the syringe does not extrude smoothly, remove the dispensing tip and extrude material directly from the syringe to make sure that material is coming out.

c. The material will start to harden if exposed to the dental light or ambient light. Be sure to protect it from light when working from a mixing pad.

d. After use, immediately remove the dispensing tip and tightly close the syringe with the wing cap.

## Clinical Hint 1

In order to inject effectively, use the surface tension of the material to restore uniformly across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

## Clinical Hint 2

Where higher abrasion resistance and thicker consistency are required, it is recommended to place GRADIA® DIRECT LoFlo instead of GRADIA® DIRECT Flo or low viscosity flowable composites.

## Clinical Hint 3

When filling a large cavity, it is recommended to place material incrementally into the cavity. Another effective method is to use GRADIA® DIRECT LoFlo for filling in undercuts or as a liner/base, and then to place composite resin (GC GRADIA® DIRECT, etc.) on top.

## 6. Light Curing

Light cure the GRADIA® DIRECT LoFlo using a light curing unit (Fig.8). Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Note :

When light curing material, wear protective glasses.

## Irradiation Time and Effective Depth of Cure

| Shade                              | A1, A2, A3 | A3.5, BW | AO3, CV |
|------------------------------------|------------|----------|---------|
| 5 sec. (Plasma arc) (2000mW/cm²)   |            |          |         |
| 10 sec. (GC G-Light) (1200mW/cm²)  | 2.0 mm     | 1.5 mm   | 1.5 mm  |
| 20 sec. (Halogen/LED) (700 mW/cm²) |            |          |         |

| Shade                              | A1, A2, A3, A3.5, AO3, CV (Cervical color), BW (Bleaching White) |
|------------------------------------|--|
| 8 sec. (Plasma arc) (2000mW/cm²)   | 3.0 mm   |
| 20 sec. (GC G-Light) (1200mW/cm²)  | 2.5 mm   |
| 40 sec. (Halogen/LED) (700 mW/cm²) | 2.0 mm   |

## 7. Shaping and Polishing

Shape and polish using standard techniques.

## SHADES

A1, A2, A3, A3.5, AO3, CV (Cervical color), BW (Bleaching White)

## STORAGE

Store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F) (Shelf life : 2 years from date of manufacture)

## PACKAGES

1. Syringe packages: syringe 1.3g (0.8mL) x 2 (in one shade), 4 dispensing tips (needle type), 2 light protective covers

2. Dispensing tip package:

- 20 dispensing tips (needle type), 2 light protective covers
- 20 dispensing tips (plastic type), 2 light protective covers

## CAUTION

1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. Take special care that the patient should not swallow material.

4. Do not mix with other similar products.

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**GRADIA® DIRECT LoFlo**  
LIGHT-CURED FLOWABLE COMPOSITE

GRADIA® DIRECT LoFlo is a multi-purpose light-cured, radiopaque, fluoride releasing, micro-filled hybrid composite resin. Thanks to its new HDR (High Density Radiopaque Pre-Polymerized Filler) technology, it contributes to achieve excellent abrasion resistance. Also the thicker consistency (lower flowability) is excellent for cavities of class III, V and root surface caries.

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**RECOMMENDED INDICATIONS**

- Restoration of Class I, II, III, IV, cavities (particularly for small Class I cavities / shallow Class V cavities / other small cavities).
- Restorations in cervical areas.
- Restorations in dentin teeth.
- Filling tunnel shaped cavities.
- Sealing hypersensitive areas.
- Uvula cements / filling in cavity undercuts.
- Sealant.
- Fixing mobile teeth.
- Additions to composite restorations. (Fig. 1)

**CONTRAINDICATIONS**

- Pulp capping.
- If any reaction to the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**DIRECTIONS FOR USE**

- 1 Preparations

1) Hold the wing cap, upright and remove the wing cap from the curing combination. Take care not to expose material to direct light from the dental lamp or natural light (Fig. 2).

2) Promptly and securely attach the dispensing tip (not the needle type) to the syringe by turning clockwise (Fig. 2). Note:

Take care not to damage the dispensing tip (too tightly) or damage its screw.

3) After attaching the dispensing tip, protect with the cover until ready for use in order to avoid exposure to light (Fig. 2).

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Note:

When using light-cured bonding system, follow manufacturer's directions for use.

3 Placement of GRADIA® DIRECT LoFlo

1) Remove cover from the dispensing tip on the syringe (Fig. 5).

2) Prior to extrusion of material, gently tap the dispensing tip outside the material to make sure that the dispensing tip is securely attached to the syringe. To remove any air from the dispensing tip, with the tip pointing upwards gently push forward until plunger reaches the mouthpiece and reaches the mouthpiece (Fig. 6).

Note:

If there is an air bubble inside the dispensing tip, air bubbles may form at the time of injection.

3) Place the dispensing tip as close as possible to the cavity, and slowly push the plunger to inject material into it (Fig. 7). Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

a. When attaching the dispensing tip, make sure that the dispensing tip is sticking to the point between the tip and the syringe in order to ensure a tight connection.

b. If the syringe does not have a stopcock, remove the dispensing tip and syringe.

Remove the dispensing tip and syringe.

c. The material will start to harden if exposed to the light.

Be sure to turn it from light when working from a mixing pad.

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Tip:

In order to effectively use the surface tension of the material to ensure separation from the dispensing tip, use a smooth surface. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface for the restoration.

**Clinical Hint 2**  
Where higher abrasion resistance and thicker consistency are required, it is recommended to place GRADIA® DIRECT LoFlo instead of GRADIA® DIRECT Flo or low viscosity flowable composites.

**Clinical Hint 3**  
When filling a large cavity, it is recommended to place material incrementally into the cavity. For example, fill one third of the cavity, then add another increment, and then place composite resin (GC GRADIA® DIRECT, etc.) on top.

**Light Curing**

Light cure the GRADIA® DIRECT LoFlo using a light curing unit (Fig. 8). Refer to the following chart for Irradiation Time and Effective Depth of Cure.

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**Irradiation Time and Effective Depth of Cure**

| Shade                             | A1, A2, A3, A5, BW | A03, CV |
|-----------------------------------|--------------------|---------|
| irradiation time                  |                    |         |
| 5 sec. (Plasma arc) (2000mW)      | 2.0 mm             | 1.5 mm  |
| 10 sec. (GC-Light) (1200mW)       |                    |         |
| 20 sec. (Halogen,LED) (700 mW/cm) |                    |         |
| 8 sec. (Plasma arc) (2000mW)      | 3.0 mm             | 2.5 mm  |
| 20 sec. (GC-Light) (1200mW)       | 3.0 mm             | 2.5 mm  |

**7. Shaping and Polishing**

Shape the GRADIA® DIRECT LoFlo using a light curing unit (Fig. 8). Refer to the following chart for Irradiation Time and Effective Depth of Cure.

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**STORAGE**

Store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F)

(Shelf life: 2 years from date of manufacture)

**PACKAGES**

1. Syringe packages: syringe 1.3g (0.8ml) x 2 (one shade), 4 dispensing tips (needle type), 2 light protective covers

2. Dispensing tip package:

• 20 dispensing tips (needle type), 2 light protective covers

• 20 dispensing tips (plastic type), 2 light protective covers

**CAUTION**

1. Material of contact with oral tissue - skin, remove immediately with cotton or a sponge and wash with alcohol. Flush with water and seek medical attention.

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5. The dispensing tip cannot be sterilized in an autoclave or chemilave.

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