




Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

1460 xxxx hotty



21-9866 18012018 / A

Replacement parts hotty, 100-120 V



Version
A

Article number
14601000



making work easy

Version	Position	Name	Article number	Wishlist
A	1	Lid HOTTY / Blue print	900034579	
A	2	Rubber bases 4 pcs.	30991	
A	101	Thermo sensing device complete	900035783	
A	102	Thermofuse complete	900035761	
A	103	Isolation Set Hotty 1457+1460	900034566	
A	104A	Printed boards 1460 (230V)	900035764	
A	104B	Printed boards 1460 (120V)	900035765	
A	105	Front label Hotty 1460	900035773	
A	106	Single-throw double-pole switch	900035747	
A	107A	Power cord (230 V)	900035762	
A	107B	Power cord (120 V)	900035763	
A	107C	Unit connector cable 230V UK	900035938	
A	107D	Power cord Argentine with cable-mount connector, straight	900032769	
A	107E	Unit connector cable 230V CH	900036070	
A	107F	Power cord Brazil with cable-mount connector, straight	900032767	
A	107G	Power cord ZA 5A with cable-mount connector, straight	900032771	
A	107H	Power cord Italy with cable-mount connector, straight	900032768	
A	108	Thermoconducting paste siliconfree	900031955	
*) Wear and tear parts, excluded from the guarantee.				

Version	Position	Name	Article number	Wishlist
A	109	Power cord 0,75mm³ w. 10A plug	900036961	
*) Wear and tear parts, excluded from the guarantee.				