




Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

5022 5000 Power pillo



21-9420 04122019 / A

Replacement parts
Power pillo



Version
A

Article number
50225000



Version	Position	Name	Article number	Wishlist
A	1A	Flat chisel, narrow, for Power pillo	50220100*	
A	1B	Flat chisel, wide, for Power pillo	50220200*	
A	1C	Channelled chisel for Power pillo	50220300*	
A	1D	Pointed chisel for Power pillo	50220400*	
A	2	Connecting hose complete	900021174	
A	3	Set of filters	900020596*	
A	4	Handpiece 2K for Power pillo	900035793	
A	5	Piston	900033307	
A	6	Power pillo one-ear clamp 11.8 – 13.8 mm	900035392	
A	101	Valve Cover	900036176	
A	102	Valve	900036175	
*) Wear and tear parts, excluded from the guarantee.				