

Renfert

making work easy



temp:ex

NEW Temporary Cement Remover

temp:ex

Temporary Cement Remover

temp:ex is a solvent for zinc oxide based, dental temporary cements.



"Until now, the removal of temporary cement has always been labor-intensive and annoying. Now, with temp:ex, this process runs simultaneously in the background. In the meantime, my assistant is able to take care of the patient."

Dr. Sigmar Schnutenhaus

Ordering information

temp:ex, 250 ml (8.45 fl.oz.)

No. 6600 0000

Recommended extension

SYMPRO, 100–240 V | 50/60 Hz

No. 6500 0000

SYMPRO mini-cup set

No. 6500 0410

Further information

www.renfert.com/temp-ex



Advantages

- Saves working time through a significant reduction in manual working steps.
- Avoids damage to the restoration via instruments during manual working steps.
- Tested compatibility with standard dental materials.

Details

- Up to 50% faster cleaning performance using temp:ex in the SYMPRO denture cleaning device.
- Economical use due to dosage aid on bottle.

