



EPOS |

An EPIC sitting experience



EPOS

When it comes to stool,
**see the benefits not
just the features.**



EPOS IS NOT JUST ANOTHER SADDLE STOOL; IT IS A TOTALLY UNIQUE AND INNOVATIVE PRODUCT WHICH:


- ✓ Designed and development based on listening to the voice of the customer to ensure that meets the needs and expectations of the user
- ✓ Made of ultra-resistant, thick, luxurious upholstery, seamless for easy & hygienic cleaning and disinfection
- ✓ Comes in 9 customizable color options
- ✓ Safe Working Load of 160kg
- ✓ Available in 2 different heights: short and medium
- ✓ 2-lever mechanism operation provides independent adjustment for the seat (height and tilt)




*Get to know the EPOS,
your personal assistant at work.*

Good posture means that the key parts of a person's body are correctly aligned and supported by the right amount of muscle tension. While working this is not always easy to remember and to maintain.


EPOS is your personal assistant that helps not only to adopt the correct sitting position, but also to maintain it all day long. How? EPOS tilts your pelvis and spine upward preventing from any abnormal position while providing at the same time a sense of comfort at your back and legs.




SIDE VIEW (INCLINATION)




TOP VIEW



SIDE VIEW



FRONT VIEW




THE EPOS STOOL IS AVAILABLE IN 2 VERSIONS:

SIZE	MINIMUM HEIGHT	MAXIMUM HEIGHT
Short (A)	50cm	62cm
Medium (A)	59cm	78cm

*** Distance of the seat surface to the floor.*

Color Shades



* Metallic Shades

FACTORY

10km. Old Highway
Thessaloniki / Thessaloniki - Kilkis
P.O. Box 142 GR 57008
T: (+30) 2310 780 692 / (+30) 2310 780 792
F: (+30) 2310 783 692
E: info@neodent.eu

ATHENS - GREECE

Feidippidou 45 (ground floor)
GR 11527, Athens
T: (+30) 210 7705020 / (+30) 2107718596
(+ 30) 210 7718595
F: (+30) 210 7718595

FOLLOW US ON



neodent
Dental equipment
www.neodent.eu